

Product data sheet Steppie

Material, both top and lower part	PPco plast – Polyprobelene copolymer
Strength	The Steppie holds 200 kilos (440 lbs).
Cleaning	Wipe off with a dry cloth or a wet one eventually added a drop of detergent
Recycling	Granulate for plast. Or incineration plant where nothing but Co ² and water are released
EU Design Registration No.	2379677
US Design Patent Registration No.	US D736,549 S
Norwegian Design Registration No.	84345
Country of origin	Denmark
Medical Device Registration	Class 1 at Danish Health and Medicines Authority, Copenhagen
Guarantee	5 years from date of delivery to the end user
Weight per piece	2.25 kilos (4.95 lbs)
Measurement per piece	560x360x53 mm (22x14x2 in)
Pieces per box	1 or 5
Weight per box	11.50 kilos (25.35 lbs)
Measurement, box, 5 pieces	365x285x565 mm (14x11x22 inches)
Measurement, box, 1 piece	365x285x60 mm (14x11x2 inches)
Boxes per pallet	4x2 boxes per layer = 8 x 3 in height = 24 boxes with 5 pieces (120 pieces per pallet)
EAN Code No., 5 pieces box	5700002000712
EAN Code No., 1 piece box	5700002000729
Steppie ApS, VAT No.	DK29823111
Customs clearance No.	9403 70 00 90, furniture of plast

Instructions of use

Each Steppie holds a removable label with reference to a QR-code and web page address with detailed instructions.

These same instructions are found on the web page (video).

Unintended effect, safety and complain

In case of any unintended effect by using the Steppie, this should immediately be reported to the dealer of the product and passed on to Steppie ApS, Denmark.

Furthermore, it is asked with regard to any kinds of safety issues or complain to pass that on to Steppie ApS as soon as possible.

Steppie ApS is committed to the *Danish Health and Medicine Authority* to file and, in certain cases, report about such incidents.